cecotec

BAMBA TOOTHCARE JET CLEAN

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SAFETY INSTRUCTIONS

Read these instructions thoroughly before using the appliance. Keep this instruction manual for future reference or new users.

- This appliance can be used by children aged 8 years and above and persons with reduced physical, sensory, or

mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved. Children must not play with the appliance. Cleaning and user maintenance should not be carried out by unsupervised children.

- This appliance should only be used with water or specific solutions.
- Make sure that the mains voltage matches the voltage stated on the appliance marking and that the wall outlet is earthed.
- Check the power cable regularly for visible damage. The power cable cannot be replaced. If damaged, discard the appliance.
- Do not operate the appliance if its cable, plug, or body show visible damage, do not work properly, or have been dropped.
- Switch off and unplug the appliance from the mains when not in use, when moving it, before assembling or removing parts, before cleaning it, and before refilling the tank. To avoid tripping hazards, wind the cable safely.
- Do not direct water under the tongue, ear, nose or other sensitive areas. Pressure generated by this appliance may cause serious damage in these areas. Read this instruction manual for proper use.
- Do not use the appliance if you are wearing oral jewellery. Remove it before use.
- Do not use the appliance if you have an open wound on your tongue or mouth. If you have been prescribed antibiotics by your physician prior to dental procedures, you should consult your dentist before using this appliance or any other oral-hygiene aid.
- Do not pull on the cable or use it to carry the appliance, do

not use the cable as a handle, do not crush the cable, and do not pull the cable around sharp edges or corners. Do not place the appliance on top of the power cable. Keep the cable away from heat sources.

- Pull from the plug, not the cable, to unplug it. Do not touch the cable or appliance with wet hands.
- Turn off the appliance before unplugging it.
- Unplug the appliance before cleaning it.
- This appliance is designed for domestic use only and is not intended for bars, restaurants, farmhouses, hotels, motels, and offices. Do not use outdoors. Do not use the appliance in places where oxygen or aerosols are being used.
- The appliance may cause bleeding gums when first used. This is completely normal and will stop happening with frequent use. Adjust the intensity of the appliance to reduce the presence of bleeding gums.
- Do not place or store the appliance where it can fall into a bathtub, shower, or sink.
- Do not immerse it in water or any other liquid. Do not use it while you shower.
- Do not use accessories not commercialised by Cecotec.
- Always follow basic safety precautions in addition to those in this manual when using electrical appliances, especially if children are present.
- It is normal for the handle to heat up during charging. However, should it heat up excessively, unplug the appliance from the mains.
- Avoid prolonged exposure of the surface of the appliance to substances such as gel, oil, soap, etc. to avoid corrosion. Disposal of the appliance and its built-in battery must comply with local regulations.
- Use the appliance only for the purposes described in this

manual. Do not use accessories that are not recommended by Cecotec.

- Never use the appliance if it has a damaged cable or plug, if it is not working properly, if it has been dropped or damaged, or if it has been dropped in water. Take the appliance to a support centre for examination and repair.
- Never block the air openings of the appliance or place it on a surface where air openings can be blocked. Keep air openings free of lint, dust, or similar.
- Do not use the appliance while sleepy.

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1. PARTS AND COMPONENTS

Fig. 1

1. Water tank
2. Electric toothbrush
3. Power/pressure-control knob
4. Suction cups
5. Tank cover/tip storage
6. Handle knob
7. Tip release button
8. Stop button
9. Main body
10. Water flosser hose

Tips. Fig. 2

- A. Regular high-pressure tip x3
- B. Periodontal tip x1
- C. Orthodontic tip x1
- D. Dental-plaque tip x1
- E. Tongue cleaner tip x1
- F. Toothbrush head x2

Note:

The graphics in this manual are schematic representations and may not exactly match the product.

2. BEFORE USE

- This appliance is packaged in a way as to protect it during transport. Take the appliance out of its box and remove all packaging materials. You can keep the original box and other packaging elements in a safe place. This will help you prevent damage to the appliance when transporting it in the future. In case the original packaging is disposed of, make sure all packaging materials are recycled accordingly.
- Make sure all parts and components are included and in good conditions. If there is any piece missing or in bad conditions, contact the official Cecotec Technical Support Service immediately.

Box content:

- Bamba ToothCare Jet Clean
- 2 toothbrush heads
- 7 water flosser tips
- Instruction manual

3. WATER FLOSSER OPERATION**Preparing the water tank**

- Remove the tank from the base. Open the cover and fill up the tank with lukewarm water. Place the tank back on the base, press firmly, and close the cover.
- The tank cover can be completely removed to make cleaning easier.

Inserting and removing tips. Fig. 3

- Insert the tip in the top-centre of the handle knob and press until you hear a click. After hearing the click, you will see the handle knob lowered and fixed in position. The coloured ring will be almost touching the end of the handle knob if the tip is correctly fixed in place.
- To remove the tip from the handle, press the tip release button and remove the tip from the handle knob.

Pressure setting

Turn the power/pressure-control knob to the lowest position for first use. Increase pressure gradually overtime until reaching the desired setting or as directed by your dental professional. Fig. 4

Using the water flosser

- Choose the tip you want to use and insert it as previously explained. Fig. 3
- Connect the plug to a power supply.
- Lean over the sink and place the tip on your mouth. Slightly close your lips to avoid splashes, but allow water to flow from your mouth to the sink.
- Point the tip at the teeth and with your free hand, turn the power/pressure-control knob clockwise until you hear a click, indicating that the appliance is in operation.
- Adjust the pressure with the power/pressure-control knob and the jet direction by turning the handle knob. Fig. 4
- Direct the water jet in a 90° angle towards the gum.
- To achieve better results, start in the molar area towards the front teeth. Fig. 5
- Slide the tip over the gumline and stop between your teeth.
- Continue the process until every area around and between teeth is clean.

Pause/Control

- You can stop water temporarily at any time by pressing the Stop button on the handle. Fig. 6
- The water flosser will automatically turn off after 3 minutes of use. Once finished, turn the appliance off by turning the power/pressure-control knob counterclockwise to the Off position and empty the tank. If you use special solutions, clean and rinse the tank thoroughly. Unplug the appliance from the mains.

Using the water flosser tips. Fig. 7**A. Periodontal tip**

- The periodontal tip releases water and anti-bacterial solution in the periodontal pockets. It is designed to reach the deepest areas in the periodontal cavity.
- To use this tip, lean over the sink and carefully place the tip under the gumline, in a 45° angle. Turn on the appliance with the lowest pressure setting and clean the gumline by putting the tip between your teeth. Continue over the gumline until you finish cleaning.

B. Regular high-pressure tip

- The high-pressure tip cleans deep inside the oral cavity with a 360-degree reach, additionally covering hard-to-clean areas, food remains, and harmful bacteria located deep inside dental crevices or under the gum.

C. Orthodontic tip

- The special orthodontic tip is exclusively designed for simultaneous brushing and rinsing of hard-to-reach areas around braces and other dental work.
- To use the orthodontic tip, lean over the sink and use the toothbrush bristles to slightly clean around the braces: water will flow through the brush for rinsing.

D. Tongue cleaner tip

- Lean over the sink, place the tongue cleaner tip on your tongue, and turn the appliance on with the lowest pressure setting. Increase it if necessary.
- Place the tip on the centre of your tongue and drag it forward while applying a slight pressure.

E. Dental-plaque tip

- Direct the jet at a 90° angle to the gum line and close the lips slightly to avoid splashing, but allow water to flow freely from the mouth into the sink. Start in the molar area (back teeth) and continue towards the front teeth, place the tip close to your teeth so bristles touch them softly.
- Gently slide the tip along the gumline, passing briefly between the teeth to brush them gently and allow water to flow between your teeth. Continue until all areas around and between the teeth have been cleaned, both inside and outside, top and bottom.

4. ELECTRIC TOOTHBRUSH OPERATION

Charging instructions

- Make sure the appliance is switched off before charging. Place the toothbrush on the charging base and make sure the latter is plugged into the mains socket. If the charging indicator light flashes, it means the toothbrush is charging. Charging is complete when the charging indicator light stops flashing and switches to steady.
- The charging indicator light may not light up instantly when the toothbrush is charged after a long period of inactivity. It should light up in 6-10 minutes.
- If the brush wasn't used for more than 6 months, charge it for at least 12 hours to reach a full charge.

Tips to extend battery life

- Completely charge the battery before first use.
- Do not recharge the brush between brushings after the first charge. Recharge it once the battery runs out.
- Completely empty the appliance at least once a year.

Using the electric toothbrush. Fig. 8

- Insert the toothbrush head onto the drive shaft.
- Wet the bristles and apply some toothpaste.
- Place the bristles against the teeth at 45 degrees to the gumline.
- Switch on the toothbrush by pressing the power/mode button, and select the mode that best suits your needs by repeatedly pressing the button.
- Move the bristles slowly in a circular motion. Do not scrub your teeth; instead, apply light pressure and let it do the brushing for you.
- Then hold down the button for 3 seconds to turn off the toothbrush. The toothbrush will automatically switch off after two minutes.
- Rinse your mouth with water after brushing.

Operating modes:

- Clean: suitable for daily toothbrushing.
- White: ideal for tooth whitening.
- Sensitive: ideal for sensitive teeth.
- Massage: to get a gentle massage while brushing.
- Gum Care: suitable for gum care.

It is recommended to replace the toothbrush head every three months.

Toothbrushing time

- The timer reminds you to change oral zones every 30 seconds for a complete brushing.
- 2 minutes is the recommended toothbrushing time.

5. CLEANING AND MAINTENANCE

Cleaning

Cleaning the water flosser

- After adding special solutions to the tank, clean it thoroughly:
- Half fill up the tank with warm tap water. With the tip pointing towards the sink, turn the appliance on and let it run until the tank is completely empty: this will help prevent the tip from clogging and help the appliance last longer.
- Unplug the appliance before cleaning the exterior surface.
- If necessary, clean the appliance with a soft cloth and a non-abrasive cleaning agent.

Cleaning the toothbrush head

Rinse the toothbrush head after every brushing and let it air dry. Do not rub the bristles to avoid deformation.

Cleaning the toothbrush handle

- Rinse and let it air dry.
- Do not immerse the toothbrush in water or other liquids.
- Do not clean the toothbrush with water above 40 °C or with abrasive cleaning agents.

Replacing the toothbrush head

- Dental professionals recommend replacing the toothbrush head every 3 months for best brushing results.
- Replace the toothbrush head as well when the bristles start to bend outwards.

Maintenance and repair

- After using the water flosser, be sure to unplug the power cable, and place the tank upside down on a flat surface. When completely dry, put it in a cool, dry place along with the body, also dry and clean. You must place it in the original box or a similar one when it is not going to be used for a long time.
- Before exposing the appliance to low temperatures, remove the tank and use the appliance until completely empty. Empty the tank also if the appliance is going to be stored for a long period of time.
- If the appliance does not work properly, contact the official Cecotec Technical Support Service.
- Do not try to disassemble or repair the appliance. If this causes any problems with the appliance, the warranty will not cover it.

6. TROUBLESHOOTING

Problem	Possible causes	Possible solutions
Leak between the tip and the handle knob.	The tip is not correctly assembled.	Remove the tip and try refitting it until it clicks.
Periodontal tip tear.	The tip is worn.	Replace it with a new tip.
The hose handle leaks.	Damaged hose.	Contact the official Cecotec Technical Support Service.
The tank leaks.	The water tank valve is missing or upside down.	Replace the tank valve. Turn the valve.
Inadequate pressure.	The tank is not assembled correctly. The tip is clogged.	Make sure the tank is full of water and correctly placed. Clean the tip following the instructions.
The unit does not release water.	The motor is not working.	Make sure the plug socket works by trying with another appliance (check the differential and/or the wall switch).
The electric toothbrush doesn't vibrate enough.	You're pressing the toothbrush too much against your teeth. Low battery level.	Do not press the bristles against your teeth. Charge the toothbrush battery.
The electric toothbrush works for a very short time.	Not charged enough. Inadequate charging.	Charge the toothbrush until the indicator light stops flashing. Make sure that the toothbrush is placed upright in the charging base and that the plug is connected.

The toothbrush does not work at all.	The toothbrush has not been used for more than 6 months. Internal mechanical failure.	Charge the toothbrush for at least 12 hours. Contact the official Cecotec Technical Support Service
--------------------------------------	--	--

7. TECHNICAL SPECIFICATIONS

Product reference: 04437
 Product: Bamba ToothCare Jet Clean
 Water flosser: 100-240 V~, 50/60 Hz, 18 W, IPX4

Electric toothbrush: Input: 5 V , 1 A, 3 W
 Battery: Li-ion 3.7 V, 800 mAh, 2.96 Wh
 IPX7

Technical specifications may change without prior notice to improve product quality.
 Made in China | Designed in Spain

8. DISPOSAL OF OLD ELECTRICAL AND ELECTRONIC APPLIANCES

This symbol indicates that, according to the applicable regulations, the product and/or batteries must be disposed of separately from household waste. When this product reaches the end of its shelf life, you should dispose of the cells/batteries/accumulators and take them to a collection point designated by the local authorities.

Consumers must contact their local authorities or retailer for information concerning the correct disposal of old machines and/or their batteries. Compliance with the above guidelines will help protecting the environment.

9. TECHNICAL SUPPORT AND WARRANTY

Cecotec shall be liable to the end user or consumer for any lack of conformity that exists at the time of delivery of the product under the terms, conditions, and deadlines established by the applicable regulations.

It is recommended that repairs are carried out by qualified personnel.

If at any moment you detect any problem with your product or have any doubt, do not hesitate to contact the official Cecotec Technical Support Service at +34 963 210 728.

10. COPYRIGHT

The intellectual property rights over the texts in this manual belong to CECOTEC INNOVACIONES, S.L. All rights reserved. The contents of this publication may not, in whole or in part, be reproduced, stored in a retrieval system, transmitted, or distributed by any means (electronic, mechanical, photocopying, recording or similar) without the prior authorization of CECOTEC INNOVACIONES, S.L.

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