cecotec

FACECARE LÜM ICE

FaceCare LÜM Ice Champagne

Máscaras LED/ LED mask

Manual de instrucciones

Instruction manual

Manuel d'instructions

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Handleiding

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SAFETY INSTRUCTIONS

Read the following instructions carefully before using the product. Keep this instruction manual for future reference or new users.

- This appliance can be used by children aged 8 years and above and people with reduced physical, sensory, or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning the use of the appliance in a safe way and understand the hazards involved. Children must not play with the appliance. Cleaning and user maintenance must not be carried out by children without supervision.
- **WARNING.** Possible eye injury (which could lead to vision loss) or skin injury if instructions are not followed. Read and follow the instructions.
- **WARNING.** Use only disposable materials and consumables recommended by the manufacturer.
- **WARNING.** Skin damage may occur after prolonged or repeated application to the same area.
- **WARNING.** Misuse can cause eye damage. Protect eyes from exposure. Do not use the appliance on eyelids or near the eyes.
- **WARNING.** Do not disable the device's built-in safety features.
- **WARNING.** If you are taking photosensitising medication or herbal remedies, consult a doctor before using the device, as it may cause skin damage.
- **WARNING.** These devices must not be used by children or on children.
- Do not use the device on areas where the integrity of the skin is compromised (e.g. in the presence of burns, injuries, blisters, scars, cuts, open wounds, active skin conditions,

recent tanning or sunburn, and/or infection).

- This device is not intended for the treatment of medical conditions. You should seek medical advice for conditions such as moles, rashes, itchy skin, fungal infections, skin growths or warts.
- Instructions on how to use it on various skin pigmentations.
- Avoid excessive use by not exceeding the specified usage parameters, such as the accumulation of pulses and/or increased frequency of use.
- Check that there is no visible damage to the emission optics or the device housing. If there is any damage, do not use the device.

Contraindications

- DO NOT use the LED Therapy Mask if you have a history of light-induced seizures. It is recommended that you test the device on a small, less sensitive area of skin (see section 5.2 – Light sensitivity test).
- DO NOT use the LED Therapy Mask to treat conditions other than those indicated. It has not been tested for other conditions and the risks are unknown.
- DO NOT use the LED Therapy Mask if you are pregnant, planning to become pregnant or breastfeeding. It has not been tested on pregnant or breastfeeding women and, therefore, the risk to the fetus, newborn or pregnant women is unknown.
- DO NOT use the LED Therapy Mask if you suffer from lupus erythematosus, photosensitive eczema, or albinism. Using the mask to treat these conditions can lead to a severe skin reaction.
- DO NOT use the LED Therapy Mask if you suffer from a photosensitivity disorder (light sensitisation). If you use

the mask and suffer from a photosensitive disorder, it can cause a severe skin reaction.

- Do not insert any objects into the device openings.
- **DONOT** use the LED Therapy Mask if you are taking medication that may cause photosensitivity. If you use the mask and are taking medication that can cause photosensitivity, it can lead to a severe skin reaction. Photosensitivity is a common side effect of various medications. These may include certain antibiotics, chemotherapy drugs, heart medications, and diuretics. If you are unsure about any medication you are taking, consult your healthcare provider.

General safety

WARNINGS:

- The LED Therapy Mask is an electrical device that, under certain circumstances, could pose a risk of electric shock to the user. Please read the attached user manual carefully and follow the indicated instructions to ensure maximum safety during its operation.
- The strap of the LED Therapy Mask contains a blend of polyester, nylon, and lycra. If you are allergic to any of these materials, you may not be able to use the strap.
- The LED Therapy Mask cannot be used when charging.
- Before each use, with the device switched off, clean the face mask.
- **DO NOT** use cleaning products other than those indicated in the cleaning instructions. The use of other cleaning products may damage the LED Therapy Mask.
- Store the LED Therapy Mask at room temperature in dry conditions, away from direct sunlight.
- **SWITCH OFF** and **REMOVE** the LED Light Therapy Mask

immediately if it feels hot or if your skin feels uncomfortably warm. Due to the length of the charging cable, strangulation may occur, so keep the product out of reach of children.

- DO NOT use the LED Light Therapy Mask without removing make-up BEFORE treatment. Some ingredients in these products can cause skin sensitisation or reduce the effect of the treatment. If you use the mask without removing these products, it may cause a skin reaction or reduce the effect of the treatment. This equipment is not intended for use by children. Children should be supervised to ensure they do not play with the equipment to avoid the risk of fire and burns.
- NEVER use the LED Light Therapy Mask near water or where water may spill onto the equipment. DO NOT use the LED Light Therapy Mask if you believe it has been in contact with water. Do not subject the equipment to strong impacts to avoid the risk of damage to the lamp.
- Keep the LED Light Therapy Mask away from HOT surfaces.
- DO NOT attempt to modify the LED Light Therapy Mask in any way. Modifying the LED Light Therapy Mask will void its warranty. Modifying the LED Light Therapy Mask could cause an electric shock when it is used again. It contains no serviceable parts. Using accessories other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity.
- DO NOT use the LED Light Therapy Mask if the mask, remote control, or cable are damaged in any way. Do not leave the equipment unattended when switched on to avoid the risk of fire or burns.
- When removing from minimum/maximum storage temperature, it is recommended to let it rest for 2 hours before use.

- Do not look directly at or direct the light into the eyes while in operation.

Electrcal safety

- DO NOT store the device at temperatures above 60°C/140°F. Increasing the internal temperature can cause malfunction or fire.
- If there is smoke or a strange odour, stop using it, switch off the device, and remove the charger.
- DO NOT damage, process, bend, twist, or heat the USB cable. Damage to the cable can cause fire or electric shock.

Meaning of symbols

Symbol	Meaning
	Refer to the instruction manual
	Class II equipment
	General warning sign
IP22	IP rating
	Keep away from sunlight
	Caution

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Remote control

Figure 2

1. LCD screen: To check the selected mode and remaining time
2. Power button
3. Mode change button: To switch between different modes.
4. Timer button: To set the usage time.

Note:

The graphics in this manual are schematic representations and may not exactly match the product.

2. BEFORE USE

- This appliance comes in a packaging designed to protect it during transport. Remove the appliance from its box. You can keep the original box and other packaging materials in a safe place to prevent damage to the appliance if you need to transport it in the future. In case the original packaging is disposed of, make sure all packaging materials are recycled accordingly.
- Check that all parts and components are included and in good condition. If any of them are missing or damaged, please contact Cecotec's Official Technical Support Service immediately.

Box contents:

- Main unit
 - Remote control
 - Type-C charging cable
 - IceTech accessory
 - This instruction manual
-
- Do not remove the product's serial number in order to keep proper traceability if technical assistance is required.

3. PRODUCT ASSEMBLY

1. Connection: Connect the mask cable to the corresponding port on the remote control.
2. Charging/Power: Ensure the device is connected to a power source or that the battery is charged before use.
3. Accessory attachment: If you wish to use the gel mask, attach it to the inside of the device using the integrated magnets.
4. Verification: Long press the power button to check that the LCD screen and LEDs activate correctly.

Note: No tools or complex assembly are required; the integrated Velcro system allows for immediate adjustment after connection.

4. OPERATION

Device setup

1. Unpack the device and all accessories.
2. Your device has an internal rechargeable battery. You will need to fully charge your device for approximately 3.5 hours before first use.

Light sensitivity test

Follow these instructions to see if you are sensitive to light:

1. Switch on the device by long pressing the "On/Off" button on the LED mask, then change colour modes with the "Mode Change" button on the remote control.
2. Place the device on a section of your forearm, maintaining contact throughout the three-minute treatment cycle.
3. Press the "colour mode selection" button again to choose a different mode for another three-minute treatment on the same forearm area.
4. Wait 24 hours. You should NOT experience any adverse reactions during the next 24 hours.
5. An adverse reaction is pain or discomfort in the treated areas or any prolonged skin redness.

Battery charging

Before use, make sure to charge the LED remote control battery as follows.

Charging time may vary depending on the battery status and can take up to 3.5 hours.

1. Connect the supplied charging cable to the Type-C port of the remote control.
2. Connect the opposite USB end to a 5 V 2 A DC charging adapter.
3. When the digital light on the remote control stops rotating and remains steady, charging is complete.

Note: The remote control indicator changes according to the charging status of the LED mask.

- **Low battery indicator:** The battery icon will flash to warn that it needs recharging. If not charged in time, the device will automatically switch off.
- **Charging indicator:** The bars of the battery icon will progressively light up on the display.
- **Full charge indicator:** The battery icon on the display will remain continuously lit.

Tips for lithium battery maintenance

- Avoid completely discharging the battery before recharging.
- Avoid overcharging the battery and limit charging time to a maximum of 8 hours.
- If the device is not used for a prolonged period, it is recommended to charge the battery at least once every three months.

How to use it

Preparation

1. Before starting, deep clean your face to remove any make-up residue or impurities. Dry your skin completely, as light technology works best on clean skin.
2. If you want a cool effect session, place the IceTech accessory, previously chilled in the refrigerator or freezer. Thanks to its magnet system, easily attach it to the inner part of the LED mask before putting it on.

Instructions for use

1. Place the mask on your face and secure it comfortably using the integrated Velcro. Ensure it stays firmly but without excessive pressure on your face.
2. Connect the remote control to the Type-C cable of the LED light therapy mask. Fig. 2
3. Press and hold the On/Off button on the remote control for 1.5 seconds. The device will automatically activate in orange light mode with a default duration of 5 minutes.
4. Press the 'M' button briefly several times in succession to cycle through the 9 light modes. Each colour has a specific function.
5. Briefly press the on/off button to switch between the 3 light intensity levels, depending on your sensitivity and needs.
6. Press the 'Time' button to select the desired session duration (5, 10, 15 or 20 minutes). The LCD screen will start a countdown to show you the time remaining.
7. The device will switch itself off once the timer reaches zero. If you wish to end the session earlier, simply press and hold the on/off button.

Tips for adjusting the headband

1. When sitting or standing, ensure that the pressure on the eye cavity is minimal, allowing the silicone pad to fit properly, while the band secures the product.
2. If the silicone pad exerts too much pressure, loosen the headband on both sides.
3. If light enters through the gap between the pad and the eye cavity, tighten the headband.

ENGLISH

4. If the mask slides, tighten the headband.
5. Select the desired mode with the buttons on the remote control.
6. If you feel discomfort, adjust the strap or remove the mask until you feel better.
7. The product will automatically switch off upon completion of the set mode and time.

Operation

The parameters for the modes are as follows:

LED Mode	Number of LEDs	Duration (minutes)	Wavelength (nm)	Light power (mW/cm ²)
ORANGE 	70 red bulbs + 70 green bulbs	5/10/15/20	630 + 530	Level 1: 10 mW/cm ² Level 2: 15 mW/cm ² Level 3: 20 mW/cm ²
CYAN 	70 blue bulbs + 70 green bulbs	5/10/15/20	465 + 530	Level 1: 15 mW/cm ² Level 2: 20 mW/cm ² Level 3: 30 mW/cm ²
WHITE 	70 red bulbs + 70 blue bulbs + 70 green bulbs	5/10/15/20	630 + 465 + 530	Level 1: 15 mW/cm ² Level 2: 20 mW/cm ² Level 3: 30 mW/cm ²
RED + NIR 	70 red bulbs + 70 infrared bulbs	5/10/15/20	630 + 850	Level 1: 10 mW/cm ² Level 2: 15 mW/cm ² Level 3: 20 mW/cm ²
NIR 	70 infrared bulbs	5/10/15/20	850	Level 1: 10 mW/cm ² Level 2: 15 mW/cm ² Level 3: 20 mW/cm ²
BLUE + NIR 	70 blue bulbs + 70 infrared bulbs	5/10/15/20	465 + 850	Level 1: 10 mW/cm ² Level 2: 20 mW/cm ² Level 3: 30 mW/cm ²

GREEN + NIR 	70 green bulbs + 70 infrared bulbs	5/10/15/20	530 + 850	Level 1: 10 mW/cm ² Level 2: 15 mW/cm ² Level 3: 25 mW/cm ²
YELLOW 	70 red bulbs + 70 green bulbs	5/10/15/20	630 + 530	Level 1: 15 mW/cm ² Level 2: 20 mW/ cm ² Level 3: 25 mW/cm ²
LILAC 	70 red bulbs + 70 blue bulbs	5/10/15/20	630 + 465	Level 1: 15 mW/cm ² Level 2: 20 mW/ cm ² Level 3: 35 mW/cm ²

Functions

- **Orange light function:** Ideal for skin with a dull or tired appearance. It helps to revitalise the face, immediately restoring natural luminosity and radiance for a "rested skin" effect.
- **Cyan light function:** It provides a soothing sensation for sensitive skin. It helps to reduce the appearance of superficial redness and helps to visually balance out excess oil.
- **White light function:** Thanks to its combination of wavelengths, it helps improve the appearance of firmness in skin lacking tone, resulting in an overall smoother, fresher and more rejuvenated look.
- **Red + NIR light function:** Our star function for mature skin. It helps to soften the appearance of expression lines and fine wrinkles, promoting skin that looks and feels more elastic and firm.
- **Infrared light function:** Boosts skin vitality. Continued use improves skin luminosity and texture, giving it a healthier and revitalised appearance.
- **Blue + NIR light function:** Designed for skin prone to imperfections. It helps keep the skin clean and looking balanced, improving the appearance of pores and preventing the appearance of oily breakouts.
- **Green + NIR light function:** The perfect product for visually evening out skin tone. It helps to minimise the appearance of dark spots and has a soothing effect that reduces signs of tiredness.
- **Yellow light function:** Focuses on improving the uniformity. It visually reduces the appearance of redness or irritation, leaving the skin feeling comforted and looking more even-toned.
- **Lilac light function:** Combines the aesthetic benefits of red and blue. It is ideal for improving the appearance of irregular skin texture and revitalising areas that have lost their natural smoothness.

Information note: This product is an aesthetic device and is not intended to diagnose, treat, cure, or prevent any disease or medical condition. It is not a medical device. If you suffer from any pre-existing skin condition, consult your dermatologist before use.

5. CLEANING AND MAINTENANCE

1. Store the LED light therapy mask at room temperature in a dry place.
2. Avoid storing it in dusty conditions, under direct sunlight, or in humid environments.
3. The LED light therapy mask contains no parts that require maintenance.
4. Before each use, with the device switched off, clean the face mask following the cleaning instructions below.

Cleaning the LED light therapy mask

1. DO NOT use other cleaning products, as they could damage the device.
2. Wipe the surface of the device with a water-based wipe or a clean, damp cloth. Stubborn stains can be gently removed with 70% alcohol. Then, use a clean, dry cloth to finish cleaning.
3. After cleaning, visually inspect the device to ensure it is clean. If it is not, repeat the cleaning steps.

Storage conditions

- Temperature: -20°C to +60°C
- Humidity: 10% to 95% RH
- Atmospheric pressure: 700 hPa to 1060 hPa

Disposal

- DO NOT dispose of any part of the LED Light Therapy Mask with household waste. The mask and remote control contain a lithium-ion polymer battery.
- To protect the environment, return this unit to the representative at the end of its useful life, or dispose of it safely following local regulations.

6. TROUBLESHOOTING

Faults and solutions

1. The remote control's LCD screen does not light up when a button is pressed, and the mask cannot be controlled.

Solution: Check if the battery isolation tab has been removed from the remote control compartment. If it has been removed and the remote still does not work, check if the battery is depleted. If it is a battery issue, you can contact the official Cecotec Technical Support Service for a replacement.

2. The remote control cannot be charged.

Solution:

- a. Replace the charging cable. If the problem is the cable, you can contact the official Cecotec Technical Support Service for a replacement.
- b. Check if the charging adapter to which the unit is connected is operational or if the voltage

indicated on the unit corresponds to the local mains voltage.

7. TECHNICAL SPECIFICATIONS

Product: FaceCare LÜM Ice / FaceCare LÜM Ice Champagne

Product reference: EU01_128883

Rated input power: 5 V $\approx$ 2 A

Rated power consumption: 7.4 W

Battery: Li-ion 3.7 V, 4000 mAh, 14.8 Wh

Maximum operating time: 3.5 hours (depending on mode)

Charging time: 3.5 hours

Protection against electric shock: Class II, Applied Part Type BF (light-emitting part)

LED mask dimensions: Approx. 69.9 x 20.8 x 0.9 (cm)

Remote control dimensions: Approx. 12.6 x 4.4 x 3.2 (cm)

LED mask weight: Approx. 250 g

Remote control weight: Approx. 140 g

Light source: Light-emitting diodes

Wavelength: Red: 630 nm; Blue: 465 nm; Green: 530 nm; Infrared: 850 nm

Light energy: Orange: 20 mW/cm²; Cyan: 30 mW/cm²; White: 30 mW/cm²; Red + NIR: 20 mW/cm²; NIR: 20 mW/cm²; Blue + NIR: 30 mW/cm²; Green + NIR: 25 mW/cm²; Yellow: 25 mW/cm²; Lilac: 35 mW/cm²

To recharge the battery, it is recommended to use an adapter with the following characteristics:

Input: 100-240 V~, 50/60 Hz

Output: 5 V $\approx$, 2 A

Technical specifications may change without prior notice to improve product quality.

Made in China | Designed in Spain

8. RECYCLING OF ELECTRICAL AND ELECTRONIC DEVICES

This symbol indicates that, according to the applicable regulations, the product and/or batteries must be disposed of separately from household waste. When this product reaches the end of its shelf life, you should dispose of the batteries/accumulators and take them to a collection point designated by the local authorities.

Consumers must contact their local authorities or retailer for information concerning the correct disposal of old electrical and electronic appliances and/or their batteries.

ENGLISH

Information regarding national packaging recycling systems and their marking can be found on our website.

Compliance with the above guidelines will help protect the environment.

9. TECHNICAL SUPPORT AND WARRANTY

Cecotec shall be liable to the end user or consumer for any lack of conformity that exists at the time of product delivery under the terms, conditions and deadlines established by applicable regulations.

Repairs should be carried out by qualified personnel.

If at any moment you detect any problem with your product or have any doubt, do not hesitate to contact the official Cecotec Technical Support Service at +34 96 321 07 28.

10. COPYRIGHT

The intellectual property rights over the texts in this manual belong to CECOTEC INNOVACIONES, S.L. All rights reserved. The content of this publication may not, either in part or in its entirety, be reproduced, stored in a retrieval system, transmitted or distributed by any means (electronic, mechanical, photocopying, recording or similar) without prior authorisation from CECOTEC INNOVACIONES, S.L.

11. SIMPLIFIED EU DECLARATION OF CONFORMITY

Cecotec Innovaciones hereby declares that this appliance complies with the essential requirements and other relevant provisions of the regulations applicable in the European Union. This appliance has been designed, manufactured and tested in compliance with required safety and quality standards. The full text of the EU Declaration of Conformity can be found on the following website: <https://cecotec.es/es/information/declaration-of-conformity>

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