cecotec

FACECARE LÜM FACECARE LÜM CHAMPAGNE

Máscaras LED/ LED mask

Manual de instrucciones
Instruction manual
Manuel d'instructions
Bedienungsanleitung
Manuale di istruzioni
Manual de instruções
Handleiding
Instrukcja obsługi
Návod k použití
Kullanma kılavuzu
Οδηγίες χρήσης
Manual d'instruccions
Használati útmutató
تامی ل عتلا ل یل

Instrucciones de seguridad	5	9. Garantie et SAV	128
Safety instructions	11	10. Copyright	128
Instructions de sécurité	16	11. Déclaration de conformité simplifiée de l'UE128	
Sicherheitshinweise	22		
Istruzioni di sicurezza	29		
Instruções de segurança	35		
Veiligheidsinstructies	41		
Instrukcje bezpieczeństwa	47		
Bezpečnostní pokyny	53		
Güvenlik talimatları	59		
Οδηγίες ασφαλείας	65		
Instruccions de seguretat	72		
Biztonsági utasítások	78		
قامال سلا تاميلا عت	84		

ÍNDICE

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SAFETY INSTRUCTIONS

Read the following instructions carefully before using the product. Keep this instruction manual for future reference or new users.

- This appliance can be used by children aged 8 years and above and people with reduced physical, sensory, or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning the use of the appliance in a safe way and understand the hazards involved. Children must not play with the appliance. Cleaning and user maintenance must not be carried out by children without supervision.
- **WARNING.** Possible eye injury (which could lead to vision loss) or skin injury if instructions are not followed. Read and follow the instructions.
- **WARNING.** Use only disposable materials and consumables recommended by the manufacturer.
- **WARNING.** Skin damage may occur after prolonged or repeated application to the same area.
- **WARNING.** Misuse can cause eye damage. Protect eyes from exposure. Do not use the appliance on eyelids or near the eyes.
- **WARNING.** Do not disable the device's built-in safety features.
- **WARNING.** If you are taking photosensitising medication or herbal remedies, consult a doctor before using the device, as it may cause skin damage.
- **WARNING.** These devices must not be used by children or on children.
- Do not use the device on areas where the integrity of the skin is compromised (e.g. in the presence of burns, injuries, blisters, scars, cuts, open wounds, active skin conditions,

- recent tanning or sunburn, and/or infection).
- This device is not intended for the treatment of medical conditions. You should seek medical advice for conditions such as moles, rashes, itchy skin, fungal infections, skin growths or warts.
 - Instructions on how to use it on various skin pigmentations.
 - Avoid excessive use by not exceeding the specified usage parameters, such as the accumulation of pulses and/or increased frequency of use.
 - Check that there is no visible damage to the emission optics or the device housing. If there is any damage, do not use the device.

Contraindications

- DO NOT use the LED Therapy Mask if you have a history of light-induced seizures. It is recommended that you test the device on a small, less sensitive area of skin (see section 5.2 – Light sensitivity test).
- DO NOT use the LED Therapy Mask to treat conditions other than those indicated. It has not been tested for other conditions and the risks are unknown.
- DO NOT use the LED Therapy Mask if you are pregnant, planning to become pregnant or breastfeeding. It has not been tested on pregnant or breastfeeding women and, therefore, the risk to the fetus, newborn or pregnant women is unknown.
- DO NOT use the LED Therapy Mask if you suffer from lupus erythematosus, photosensitive eczema, or albinism. Using the mask to treat these conditions can lead to a severe skin reaction.
- DO NOT use the LED Therapy Mask if you suffer from a photosensitivity disorder (light sensitisation). If you use

the mask and suffer from a photosensitive disorder, it can cause a severe skin reaction.

- Do not insert any objects into the device openings.
- **DONOT** use the LED Therapy Mask if you are taking medication that may cause photosensitivity. If you use the mask and are taking medication that can cause photosensitivity, it can lead to a severe skin reaction. Photosensitivity is a common side effect of various medications. These may include certain antibiotics, chemotherapy drugs, heart medications, and diuretics. If you are unsure about any medication you are taking, consult your healthcare provider.

General safety

WARNINGS:

- The LED Therapy Mask is an electrical device that, under certain circumstances, could pose a risk of electric shock to the user. Please read the attached user manual carefully and follow the indicated instructions to ensure maximum safety during its operation.
- The strap of the LED Therapy Mask contains a blend of polyester, nylon, and lycra. If you are allergic to any of these materials, you may not be able to use the strap.
- The LED Therapy Mask cannot be used when charging.
- Before each use, with the device switched off, clean the face mask.
- **DO NOT** use cleaning products other than those indicated in the cleaning instructions. The use of other cleaning products may damage the LED Therapy Mask.
- Store the LED Therapy Mask at room temperature in dry conditions, away from direct sunlight.
- **SWITCH OFF** and **REMOVE** the LED Light Therapy Mask

immediately if it feels hot or if your skin feels uncomfortably warm. Due to the length of the charging cable, strangulation may occur, so keep the product out of reach of children.

- DO NOT use the LED Light Therapy Mask without removing make-up BEFORE treatment. Some ingredients in these products can cause skin sensitisation or reduce the effect of the treatment. If you use the mask without removing these products, it may cause a skin reaction or reduce the effect of the treatment. This equipment is not intended for use by children. Children should be supervised to ensure they do not play with the equipment to avoid the risk of fire and burns.
- NEVER use the LED Light Therapy Mask near water or where water may spill onto the equipment. DO NOT use the LED Light Therapy Mask if you believe it has been in contact with water. Do not subject the equipment to strong impacts to avoid the risk of damage to the lamp.
- Keep the LED Light Therapy Mask away from HOT surfaces.
- DO NOT attempt to modify the LED Light Therapy Mask in any way. Modifying the LED Light Therapy Mask will void its warranty. Modifying the LED Light Therapy Mask could cause an electric shock when it is used again. It contains no serviceable parts. Using accessories other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity.
- DO NOT use the LED Light Therapy Mask if the mask, remote control, or cable are damaged in any way. Do not leave the equipment unattended when switched on to avoid the risk of fire or burns.
- When removing from minimum/maximum storage temperature, it is recommended to let it rest for 2 hours before use.

- Do not look directly at or direct the light into the eyes while in operation.

Electrical safety

- DO NOT store the device at temperatures above 60°C/140°F. Increasing the internal temperature can cause malfunction or fire.
- If there is smoke or a strange odour, stop using it, switch off the device, and remove the charger.
- DO NOT damage, process, bend, twist, or heat the USB cable. Damage to the cable can cause fire or electric shock.

Meaning of symbols

Symbol	Meaning
	Refer to the instruction manual
	Class II equipment
	General warning sign
IP22	IP rating
	Keep away from sunlight
	Caution

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1. PARTS AND COMPONENTS

Figure 1

1. Face mask with LED technology
2. Headband
3. Remote control
4. USB cable
5. Charging cable connection ports

Note:

The graphics in this manual are schematic representations and may not exactly match the product.

2. BEFORE USE

- This device comes in packaging designed to ensure it remains protected during transport. Remove the product from its box and, if possible, keep the original packaging in a safe place for future transport. If you decide to dispose of it, please do so responsibly by following local recycling regulations.
- Check that all components are included and in perfect condition. If you notice any faults or find that any accessories are missing, please contact the official Cecotec Technical Support Service immediately.

Box contents:

- Flexible silicone face mask + connecting cable
 - Headband (x2)
 - Compact remote control
 - USB charging cable
 - This instruction manual
-
- For warranty and traceability purposes, do not remove the serial number affixed to the device. This will be essential when requesting technical support.

3. OPERATIONAL SAFETY

The LED technology mask is a precision electronic device. To ensure a safe experience and avoid the risk of electric shock, it is essential to follow these guidelines:

- **Must read:** Read this instruction manual carefully before your first session and strictly follow the operating instructions.
- **Health:** The user must be in good general health to use the device.
- **Safety warning:** As with any electrical equipment, misuse or use in unsuitable circumstances (such as contact with water) could pose a risk to the user.

4. OPERATION

1. **Preparation:** Before you begin, ensure that the remote control battery is fully charged. Cleanse your face thoroughly to remove any traces of make-up or impurities; this will help the light to penetrate more effectively.
2. **Connection:** Connect the mask to the remote control using the supplied USB cable.
3. **Placement:** Place the mask on your face and adjust it comfortably using the elastic headbands.
4. **Activation:** Press and hold the power button to start the device.
5. **Mode selection:** Use the remote control to choose between Automatic Mode or Manual Mode (selecting the desired light colour according to your skin type and aesthetic preferences).
6. **Session:** Enjoy the beauty treatment for up to 15 minutes. The mask will switch off automatically once the time has elapsed.
7. **Finishing:** Remove the mask and proceed with your usual moisturising or serum routine.

Intended use of the 9 lights

- **Red light + IR (630 nm + 830 nm):** This is the go-to combination for comprehensive anti-ageing care. Red light helps the skin appear firmer, whilst near-infrared (NIR) light penetrates to boost vitality. It is ideal for improving the elasticity and softening fine lines following your beauty routine.
- **Blue light + IR (415 nm + 830 nm):** Specifically designed to balance skin prone to blemishes. Whilst blue light helps keep the skin's surface clean and free from impurities, infrared soothes the feeling of tightness. It is recommended for oily skin seeking a more matte and uniform appearance.
- **Green light + IR (520 nm + 830 nm):** This combination specialises in brightening the face and improving the appearance of uneven pigmentation. Green light works to promote a more even tone, making it perfect for achieving radiant skin and visually minimising the appearance of sunspots.
- **Near-infrared light or NIR (830 nm):** This light is invisible to the human eye, but acts on the surface layers to provide deep revitalisation. Its main benefit is to improve the sensation of elasticity and enhance a healthy, rested appearance from within.
- **Lilac light (630 nm + 415 nm):** Acts as a synergistic support for the appearance of the skin barrier. It helps improve skin texture and provides soothing care, making it ideal for skin showing signs of old blemishes or sensitive areas.

ENGLISH

- **Yellow light (630 nm + 520 nm):** Its main function is to visibly reduce redness and enhance the skin's comfort. It is the preferred choice for easily irritated skin, helping the face to appear more balanced and less congested.
- **White light (630 nm + 415 nm + 520 nm):** This is the 'all-in-one' daily care spectrum. By combining the benefits of the primary lights, it offers a comprehensive experience that addresses radiance, firmness and clarity in a single step, ideal for general facial beauty maintenance.
- **Orange light (630 nm + 520 nm):** Focuses on revitalising a dull complexion. It is excellent for improving the appearance of tired-looking skin (such as after a poor night's sleep), whilst also helping facial care products (creams and serums) to be better absorbed.
- **Cyan light (415 nm + 520 nm):** Provides gentle care for skin prone to shine, with a highly relaxing effect. By combining blue light with the soothing benefits of green light, it is the preferred choice for improving the appearance of enlarged pores on sensitive skin.

Information note: This product is an aesthetic device and is not intended to diagnose, treat, cure, or prevent any disease or medical condition. It is not a medical device. If you suffer from any pre-existing skin condition, consult your dermatologist before use.

Product features

- **Premium, minimalist design:** High-end aesthetics designed for superior ergonomic fit to the face, optimising the reach of the LED light.
- **Advanced eye protection:** Includes high-quality silicone pads that protect the eyes from the intensity of the light.
- **Fully customisable:** Allows you to select the mode of use and adjust the duration according to your preferences and skin goals.
- **Portability and comfort:** Compact and lightweight design for easy transport, enabling advanced facial care anytime, anywhere.

Remote control functions

The remote control manages the on/off function and powers the mask via its rechargeable lithium-ion polymer battery. The device has two operating modes:

- Automatic mode (A): Allows you to enjoy a full 15-minute session. This mode automatically alternates between red, blue and green lights at 5-minute intervals, keeping the infrared light active throughout the process for a complete revitalisation of the skin's appearance.
- Manual mode (M): Allows you to select specific colours according to your desired beauty goal. After each selection, the device will operate for a maximum of 15 minutes.

Mode and colour guide

Indicator	Treatment description
R + I	Red mask light + NIR
B + I	Blue mask light + NIR
G + I	Green mask light + NIR
NIR	Near-infrared or NIR mask light
R + B	Lilac mask light
R + G	Yellow mask light
R + B + G	White mask light
R + G	Orange mask light
B + G	Cyan mask light
A	Automatic mask operating mode
M	Manual mask operating mode

Contraindications, warnings and precautionsContraindications

- DO NOT use the mask to treat conditions other than those indicated. It has not been tested for other conditions and the risks are unknown.
- DO NOT use the mask if you are pregnant, planning to become pregnant or breastfeeding. It has not been tested in these circumstances, so the risk is unknown.
- DO NOT use the mask if you suffer from lupus erythematosus, photosensitive eczema or albinism, as it could cause a severe skin reaction.
- DO NOT use the mask if you suffer from any photosensitive disorder. It could cause a severe skin reaction.
- DO NOT use the mask if you are taking medicines that cause photosensitivity, such as certain antibiotics, chemotherapy drugs and diuretics. Consult your doctor if you have any doubts.
- Avoid using it on wounds.
- The user is the intended operator and can safely use the skincare function.
- Do not carry out maintenance whilst the product is in use, including charging.
- Use the mask with caution.

Warnings:

- Carry out a sensitivity test before using the mask. DO NOT use it if you are sensitive to it.

ENGLISH

- DO NOT use the mask in any way other than that specified in the instructions. Untested use may be dangerous.
- SWITCH OFF and REMOVE the mask if it feels uncomfortably warm.
- Use only the recommended hydrogel mask. Other masks have not been tested.
- Remove creams, lotions or make-up BEFORE using the mask to avoid sensitisation or reduced effectiveness.
- DO NOT use the mask near water or if it has been in contact with water.
- Keep the mask away from hot surfaces.
- DO NOT attempt to modify the mask, as this will void the warranty and may cause an electric shock.
- DO NOT use the mask if it is damaged.
- DO NOT use it whilst charging, and do not touch the device with wet hands whilst charging.
- DO NOT use the product in areas where medicines are stored.
- Avoid storing flammable or explosive materials near the product.
- Store the product in a cool, dry place; do not refrigerate it.
- DO NOT use the product in bathrooms, swimming pools or other damp places.
- Keep the product away from strong magnetic fields such as those from microwaves or televisions.
- DO NOT place the product near other equipment without checking that it is working properly.
- Keep the product out of the reach of children to prevent misuse.
- DO NOT expose the product to the open air for prolonged periods.
- Avoid bending or applying pressure to the product and the power cable.
- DO NOT modify or replace the power cable and battery yourself.
- DO NOT look directly at the light source of the LED mask.
- Avoid prolonged use in areas with unstable voltage.
- DO NOT store magnets or metals near the product.
- Avoid exposure to harmful optical radiation.
- DO NOT touch cables and hoses to avoid strangulation.
- Avoid inhaling or swallowing small parts.
- People allergic to the materials should not use the equipment.
- DO NOT touch the device with wet hands to avoid injury.
- Use only accessories and materials described in the instructions.
- DO NOT connect this device to others not specified in the instructions.

Precautions

- Use only cables and accessories sold by our company to avoid increasing the LED therapy mask's emission.
- Install and use the product in accordance with the instructions in the manual.
- This product is for personal use only; sharing it is not recommended.

- The headband contains polyester and nylon. If you are allergic, avoid using it.
- The user must be trained in the proper use of the device.
- The remote control is Class II for protection against electric shock; the mask has internal protection.
- Report any risks of interference during specific treatments.
- If connected to a separate power supply, this must be specified as part of the system.
- Protect your eyes and skin; do not look directly at the light.
- The skin care function is safe; avoid using unrelated functions.
- Do not carry out maintenance yourself, such as changing batteries.
- Using unspecified controls or settings may result in exposure to hazardous radiation.

Operating mode information

This device uses precision light technology that combines three core wavelengths to create nine specific beauty programmes. Each mode is designed to improve the appearance of the skin through the controlled emission of light:

- **Red light (630 nm):** Anti-ageing action and improvement in the appearance of firmness.
- **Blue light (415 nm):** Helps keep the skin clear and balanced, reducing the appearance of blemishes.
- **Green light (520 nm):** Helps to visually even out skin tone and adds radiance.
- **Near-infrared light or NIR (830 nm):** Deep revitalisation for healthier-looking skin (light invisible to the human eye).

Combined programmes (light synergy)

The mask optimises results by simultaneously combining light spectra:

- **Rejuvenation mode (red + NIR):** Ideal for softening the appearance of fine lines and improving the feeling of smoothness.
- **Purifying mode (blue + NIR):** Pore care and a soothing sensation after the appearance of shine or blemishes.
- **Brightening mode (green + NIR):** Improves the appearance of superficial blemishes and revitalises the face.
- **Recovery mode (lilac):** A fusion of red + blue to visually improve the texture of skin with old marks.
- **Anti-redness mode (yellow):** A fusion of red + green to soothe the skin and enhance the feeling of skin protection.
- **Vitality mode (orange):** A fusion of red + yellow to revitalise the skin and enhance its radiant appearance.
- **Balance mode (cyan):** A combination of blue and green for a soothing effect on skin prone to sensitivity.
- **Total mode (white):** Triple spectrum (red + blue + green) for comprehensive maintenance of the face's beauty.

ENGLISH

Sensitivity test

Before using the device for the first time, carry out a sensitivity test by following these steps:

1. Place the device on the inside of your arm.
2. Switch on the remote control and let it run for the full 10-minute cycle.
3. If at any point you feel heat or discomfort, SWITCH OFF the device immediately.
4. Wait at least 6 hours after treatment. If your skin becomes red or itchy, the device is not suitable for you. Return it to the retailer within 14 days of purchase.

How to use the LED light therapy mask

Recommendations

For best results, carry out sessions 1–2 times a week for a minimum of 4 weeks. The maintenance routine should be carried out as needed, depending on your lifestyle and skin condition.

Performing a beauty session

1. Assemble your mask. Attach the headbands to the mask, ensuring they are securely fastened. Fig. 2

Note: Ensure the smooth sides are in contact with the skin and the textured sides face outwards for a better fit.

Figure 3

1. Soft side
 2. Rough side
 3. Join both sides
2. Clean and dry your face and hands thoroughly.
 3. Fitting the mask: Ensure that the inner side of the mask is in contact with the skin and adjust the headbands for a comfortable fit.

Figure 4

- A. Outer side
 - B. Inner side
4. Adjust the mask's headband to ensure a comfortable fit.
 5. Sit or lie down and place the mask over your face.
 6. Use the remote control to select the light colour, brightness, flashing or steady light, and session duration.
 7. Apply your serums or skincare creams after the treatment, as their absorption is enhanced following the use of light therapy.
 8. Continue with your skincare routine. Carry out 15-minute sessions, 3–4 times a week, and

reduce to 1–2 times a week once the desired results are achieved.

Tips

- In manual mode, ensure the mask is used for no longer than 15 minutes to avoid skin damage.
- If this is your first time using the product, follow this usage schedule and adjust it according to your skin condition:
 - First week: Use once every 3 days, for 10 minutes per session.
 - Second week: Use once every 2 days, 10 minutes per session.
 - Third week: Use once every 2 days, 10-15 minutes per session.
 - Fourth week: Use daily, 10-15 minutes per session.

Using the remote control

Modes

The mask offers manual and automatic modes.

Manual mode

- Press and hold the button on the mask for one second to switch it on; it will enter manual mode by default (red light).
- If the indicator light flashes 5 times without any operation, the mask will switch off automatically after 10 minutes.

Automatic mode

- Press the mode button once to switch to automatic mode.
- The device will start a 15-minute comprehensive care cycle, designed to provide a complete and balanced experience for your skin without the need for manual adjustments. The mask will follow this automatic sequence:
 1. 5 minutes: Red light + NIR (rejuvenation mode).
 2. 5 minutes: Blue light + NIR (purifying mode).
 3. 5 minutes: Green light + NIR (brightening mode).
- The mask will switch off automatically at the end of the 15-minute cycle. You can also switch it off manually by holding down the power button for 3 seconds.

Light selection in manual mode

1. Press and hold the power button for 1 second to switch on the mask. The red light and the indicator light on the left-hand side of the button will light up. To switch it off, press and hold the button for 3 seconds or press it four times.
2. Press the power button again to activate the blue light. The indicator light at the top left of the button will light up. To switch it off, press and hold the button for 3 seconds or press it three times.

ENGLISH

3. Press the power button a third time to activate the infrared light. The indicator light at the top right of the button will light up. To switch it off, press and hold the button for 3 seconds or press it two times.
4. Press the power button a fourth time to activate the mixed light. The indicator light on the right-hand side of the button will light up. To switch it off, press and hold the button for 3 seconds or press it once.

Light selection in automatic mode

1. With the mask switched on, press the auto/manual button to switch between manual (M) and automatic (A) mode.
2. In automatic mode, if there is no activity for 3 seconds, the timer will start.
3. It will switch off automatically after 15 minutes.

Figure 5 shows how to use the remote control.

Charging

Fig. 6

Important: Before charging, make sure you have the mask, the USB cable and the remote control ready.

1. The mask comes with a remote control.
2. Connect the USB cable to the charging port on the remote control.
3. Plug the USB adapter into the other end, ensuring that the product's voltage matches that of the power supply.
4. During charging, the indicator on the remote control will light up white. All indicators will light up when charging is complete.

Note: Battery type: lithium-ion

Warnings:

- Avoid using this equipment near other devices or stacked with them, as this may affect its operation. If necessary, monitor both pieces of equipment to ensure normal operation.
- The use of accessories, transducers and cables not specified by the manufacturer may increase electromagnetic emissions or reduce the equipment's electromagnetic immunity, affecting its operation.
- Do not use portable RF communications equipment within 30 cm of any part of the equipment, including cables specified by the manufacturer, to avoid degrading its performance.

Tips:

- When charging, use only the original USB cable supplied with the mask and do not use the mask whilst it is charging.
- As this product has a battery capacity of 1800 mAh, the adapter is not included in the standard accessories. We recommend charging with a 0.5 A adapter to reduce charging time.
- When the battery voltage is below 3.2 V, the indicator light on the remote control will flash and the mask will switch off automatically after flashing 10 times.

Risks and benefitsKnown risks

- Immediately after beauty sessions, the skin may appear red. The skin should return to normal within 1-2 hours.
- Temporary headaches immediately following light treatment. Headaches can be treated with a suitable painkiller and rehydration.
- Reduces fine lines on the face and improves skin texture and pigmentation.
- Reduces skin imperfections.

Environmental conditions

Item	Specifications	
Operating environment	Temperature	5°C ~ +35°C
	Humidity	85% RH or less
	Atmospheric pressure	80-106 kPa
	Altitude	≤ 2000 metres
Storage environment	Temperature	-10°C ~ +50°C
	Humidity	85% RH or less
	Atmospheric pressure	80-106 kPa
	Altitude	≤ 2000 metres

Environmental conditions the product can withstand

- The mask can be placed in an environment with a temperature of $60 \pm 2^\circ\text{C}$ and a humidity of $85\% \pm 2\%$ RH for 2 hours, and then in an environment with a temperature of $-20 \pm 2^\circ\text{C}$. After repeating the above 9 times, the mask is placed at room temperature for 2 hours.

ENGLISH

- The mask can be placed in an environment with a temperature of $60 \pm 2^{\circ}\text{C}$ and a humidity of $65\% \pm 2\% \text{ RH}$ for 48 hours, and then at room temperature for 2 hours.
- The mask may be placed in an environment with a temperature of $-20 \pm 2^{\circ}\text{C}$ for 48 hours, and then at room temperature for 2 hours.
- The mask may be placed in an environment with a temperature of $40 \pm 3^{\circ}\text{C}$ and a humidity of $90\% \pm 2\% \text{ RH}$ for 48 hours, and then at room temperature for 2 hours.
- The mask may be placed in an environment with a temperature of 40°C for 2 hours, and the relevant functions tested in this environment.
- The mask may be placed in an environment with a temperature of 0°C for 2 hours, and the relevant functions tested in this environment.

5. CLEANING AND MAINTENANCE

- Clean the LED mask after each use.
- Store it at room temperature in a dry place.
- Avoid dusty areas, direct sunlight and humidity.
- The mask requires no internal maintenance.
- Before each use, clean the mask following these instructions.

Cleaning the LED mask

6. Switch off the device before cleaning.
7. Use a 70%-80% ethanol solution to wipe the surface twice for 3 minutes, then wipe with a damp or dry cloth.
8. Do not use benzene, thinners, butyl acetate or other flammable materials.
9. Do not spray water directly onto the device.
10. Avoid other cleaning products that may damage the plastics.

6. TROUBLESHOOTING

Problem	Possible cause	Action
The device does not switch on.	1. The plug is not properly inserted. 2. The LED mask is not properly connected. 3. The power button was not presses and held down for one second.	1. Insert the plug correctly. 2. Check the power cable. 3. Press and hold the power button again.

<p>The LED light does not turn on.</p>	<ol style="list-style-type: none"> 1. The LED mask is not properly connected. 2. The circuit board is faulty. 3. The power cable has been pulled too hard. 	<ol style="list-style-type: none"> 1. Check the power cable. 2. Do not pull too hard on the power cable.
<p>Flashing during use. The light mode cannot be changed.</p>	<ol style="list-style-type: none"> 1. The remote control's battery voltage is too low. 2. The power cable is not fully inserted into the remote control's port. 3. The device is faulty. 	<ol style="list-style-type: none"> 1. Charge the remote control. 2. Check that the power cable is in good condition. 3. Contact the official Cecotec Technical Support Service.

- If the problem persists, contact the official Cecotec Technical Support Service.
- DO NOT attempt to repair the appliance yourself.
- DO NOT modify the device. The LED system has NO user-serviceable parts.

7. TECHNICAL SPECIFICATIONS

Product: FaceCare LÜM / FaceCare LÜM Champagne

Product reference: EU01_124613

Voltage: 5 V =

Intensity: 1 A

Power: 5 W

Battery: Li-ion 3.7 V, 1800 mAh, 6.66 Wh

Gross weight: 250 g

LED mask dimensions: approx. 207×277×43 mm

To recharge the battery, it is recommended to use an adapter with the following characteristics:

Input: 100-240 V~, 50/60 Hz

Output: 5 V =, 1 A

Technical specifications may change without prior notice to improve product quality.

Made in China | Designed in Spain

8. RECYCLING OF ELECTRICAL AND ELECTRONIC DEVICES

This symbol indicates that, according to the applicable regulations, the product and/or batteries must be disposed of separately from household waste. When this product reaches the end of its shelf life, you should dispose of the batteries/accumulators and take them to a collection point designated by the local authorities.

Consumers must contact their local authorities or retailer for information concerning the correct disposal of old electrical and electronic appliances and/or their batteries.

Information regarding national packaging recycling systems and their marking can be found on our website.

Compliance with the above guidelines will help protect the environment.

9. TECHNICAL SUPPORT AND WARRANTY

Cecotec shall be liable to the end user or consumer for any lack of conformity that exists at the time of product delivery under the terms, conditions and deadlines established by applicable regulations.

Repairs should be carried out by qualified personnel.

If at any moment you detect any problem with your product or have any doubt, do not hesitate to contact the official Cecotec Technical Support Service at +34 96 321 07 28.

10. COPYRIGHT

The intellectual property rights over the texts in this manual belong to CECOTEC INNOVACIONES, S.L. All rights reserved. The content of this publication may not, either in part or in its entirety, be reproduced, stored in a retrieval system, transmitted or distributed by any means (electronic, mechanical, photocopying, recording or similar) without prior authorisation from CECOTEC INNOVACIONES, S.L.

11. SIMPLIFIED EU DECLARATION OF CONFORMITY

Cecotec Innovaciones hereby declares that this appliance complies with the essential requirements and other relevant provisions of the regulations applicable in the European Union. This appliance has been designed, manufactured and tested in compliance with required safety and quality standards. The full text of the EU Declaration of Conformity can be found on the following website: <https://cecotec.es/es/information/declaration-of-conformity>

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