cecotec

HEATCONFORT THERMALNECK

Almohadilla eléctrica cuello/ Electric heating pad for the neck

Manual de instrucciones
Instruction manual
Manuel d'instructions
Bedienungsanleitung
Manuale di istruzioni
Manual de instruções
Handleiding
Instrukcja obsługi
Návod k použití
Kullanma kılavuzu
Οδηγίες χρήσης
Manual d'instruccions
Használati útmutató
دليل التعليمات

Instrucciones de seguridad	5
Safety instructions	8
Instructions de sécurité	12
Sicherheitshinweise	16
Istruzioni di sicurezza	20
Instruções de segurança	23
Veiligheidsinstructies	27
Instrukcje bezpieczeństwa	30
Bezpečnostní pokyny	34
Güvenlik talimatları	37
Οδηγίες ασφαλείας	41
Instruccions de seguretat	45
Biztonsági utasítások	48
52	قَدَمُ السَّلَامَاتِ

ÍNDICE

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SAFETY INSTRUCTIONS

Read the following instructions carefully before using the product. Keep this instruction manual for future reference or new users.

- This appliance is intended for domestic use only and is not suitable for use in commercial establishments such as bars, restaurants, farms, hotels, motels, and offices.
- If the power cord is damaged, it must be replaced by the manufacturer, its after-sales service, or similarly qualified personnel to avoid a hazard.
- Children under 3 years of age should not use the appliance, as they are not capable of reacting to overheating.
- Children aged 3 to 8 may only use the appliance under supervision and must always set the temperature to the minimum setting.
- This appliance can be used by children aged 8 years and above and people with reduced physical, sensory, or mental capabilities or lack of experience and knowledge if

they have been given appropriate supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved. Children must not play with the appliance. Cleaning and user maintenance must not be carried out by children without supervision.

- Do not use the appliance on children or people with reduced sensitivity to heat or who are asleep, as they do not perceive heat in the same way.
- Before each use, check the heating pad for any signs of damage or defects. The appliance must not be used if there are signs of damage.
- The appliance must not be used for heating animals.
- Do not operate the heating pad if the controller or the power cord are damaged or if you notice that the appliance is not working properly.
- Before connecting the appliance to the mains, make sure that the voltage indicated on the technical specification label corresponds to the voltage in your home.
- Do not bend or crumple the heating pad during use.
- Do not stick needles, pins, safety pins or other sharp objects into the appliance.
- Avoid falling asleep while the heating pad is on.
- Prolonged use of the appliance at maximum temperature may cause skin burns.
- When storing, do not place any objects on top that could cause creases; allow to cool before folding. The controller should never be covered or placed on or under the heating pad while it is in operation.
- Do not sit, lie down or put your body weight on the heating pad. Place it on the back of your neck and spread the flaps over your shoulders, over a thin, smooth garment, ensuring even contact. It is designed exclusively for the neck and shoulders.

- Do not touch the appliance if it has fallen into water. Unplug it immediately from the socket.
- Keep the cable away from hot surfaces.
- This device is not intended for medical use in hospitals.
- Do not pull on the cable or the appliance to disconnect the plug from the socket.
- Do not use the heating pad when it is wet or damp. Only use it in dry rooms (avoid using it in the bathroom or similar places).
- Do not expose the controller and cable to weather or moisture.
- In case of damage, do not try to repair the appliance on your own. Repairs must only be carried out by qualified personnel.
- The heating pad should only be used with the controller provided with the appliance.
- If you have any health problems, consult your doctor before using the heating pad.
- Do not use the heating pad on parts of the body that are inflamed, injured or swollen.
- If you experience constant pain in muscles and joints, inform your doctor. This type of pain can be a symptom of a serious illness.
- If you experience any pain or discomfort while using the appliance, stop using it immediately.

Symbol	Meaning
	Read the instruction manual.
	Class II appliance

	Do not stick needles or pins into the device
	Do not use if bent or folded
	Not suitable for children under 3
	For domestic use only
	Normal wash, max. 30°C
	Do not use bleach
	Do not tumble dry
	Do not iron
	Do not dry clean

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Note:

The graphics in this manual are schematic representations and may not exactly match the product.

2. BEFORE USE

- This appliance comes in a packaging designed to protect it during transport. Remove the appliance from its box. You can keep the original box and other packaging materials in a safe place to prevent damage to the appliance if you need to transport it in the future. In case the original packaging is disposed of, make sure all packaging materials are recycled accordingly.
- Check that all parts and components are included and in good condition. If any of them are missing or damaged, please contact the official Cecotec Technical Support Service immediately.

Box contents:

- Product
- Instruction manual
- Do not remove the product's serial number in order to keep proper traceability if technical assistance is required.

3. OPERATION

- The neck heating pad provides localised heat to the nape of the neck, the trapezius muscles and the upper back. Place it so that the top edge rests at the nape of your neck and spread the flaps over your shoulders. Adjust the front fastening to keep it in place. Use it on a thin garment and ensure that the surface is completely flat, without any creases, and in even contact. Keep the throat area clear.
- Temperature setting: use the temperature setting button to choose from 10 settings to achieve your desired level of comfort within the range (30°C to 65°C). The current setting is shown on the digital display of the controller.

ENGLISH

- Timer: Use the Timer button to select a running time between 1 and 6 hours. This function promotes energy saving and helps reduce the risk of skin overheating.

Instructions for use

1. First connect the power cable to the pad and then to a mains socket. The pad enters standby mode.
2. Press the on/off button to start the operating mode.
3. Select the temperature: press the temperature setting button to choose the temperature.
4. Activate the timer: press the timer button to select a running time between 1 and 6 hours.

Important

- To prevent **error code F2** from appearing, first **connect the cable to the device** and **then to the mains socket**.
If **error F2** still appears, **disconnect the product from the mains**, wait a few seconds and **reconnect it** in the same order (first to the device and then to the mains).
- When the appliance is switched on, the default settings are: 6-hour timer (auto-off) and temperature set to Level 1 (30°C).
- Use the controller only with this neck heating pad. Do not connect it to any other heating devices.

4. CLEANING AND MAINTENANCE

Cleaning

1. Switch the unit off using the power button.
 2. Unplug the power cord from the mains socket.
 3. Remove the controller and power supply unit.
- It is recommended to wash the neck heating pad by hand and reassemble it only when it is completely dry.
 - Do not use chemical cleaners.
 - Air-dry the pad by laying it flat on a well-ventilated surface, away from direct sunlight and heat sources (heaters, stoves, tumble dryers).
 - Do not iron, wring out or squeeze the neck heating pad. Do not tumble dry.
 - Do not reconnect the power cable until the pad is completely dry.

Storage

After use, switch off the pad, unplug the power cable from the mains socket and store the product in a dry, safe, clean and cool place.

5. TECHNICAL SPECIFICATIONS

Product: HeatConfort ThermalNeck

Product reference: EU01_133872

Voltage: 230 V~

Frequency: 50 / 60 Hz

Power: 100 W

Technical specifications may change without prior notice to improve product quality.

Made in China | Designed in Spain

6. RECYCLING OF ELECTRICAL AND ELECTRONIC DEVICES

This symbol indicates that, according to the applicable regulations, the product and/or batteries must be disposed of separately from household waste. When this product reaches the end of its shelf life, you should dispose of the batteries/accumulators and take them to a collection point designated by the local authorities.

For detailed information on how to properly dispose of electrical and electronic appliances and/or batteries, consumers should contact their local authorities.

Information regarding national packaging recycling systems and their marking can be found on our website.

Compliance with the above guidelines will help protect the environment.

7. TECHNICAL SUPPORT AND WARRANTY

Cecotec shall be liable to the end user or consumer for any lack of conformity that exists at the time of product delivery under the terms, conditions and deadlines established by applicable regulations.

Repairs should be carried out by qualified personnel.

If you detect any problem with the product or have any queries, please contact the official Cecotec Technical Support Service at +34 96 321 07 28.

8. COPYRIGHT

The intellectual property rights over the texts in this manual belong to CECOTEC INNOVACIONES, S.L. All rights reserved. The content of this publication may not, either in part or in its entirety, be reproduced, stored in a retrieval system, transmitted or distributed by any means (electronic, mechanical, photocopying, recording or similar) without prior authorisation from CECOTEC INNOVACIONES, S.L.

9. SIMPLIFIED EU DECLARATION OF CONFORMITY

 Cecotec Innovaciones hereby declares that this appliance complies with the essential requirements and other relevant provisions of the regulations applicable in the European Union. This appliance has been designed, manufactured and tested in compliance with required safety and quality standards. The full text of the EU Declaration of Conformity can be found on the following website: <https://cecotec.es/es/information/declaration-of-conformity>

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