cecotec

SKINCARE MASSAGE GUN

Pistola masajeadora / Massage gun

Manual de instrucciones
Instruction manual
Manuel d'instructions
Bedienungsanleitung
Manuale di istruzioni
Manual de instruções
Instructiehandleiding
Instrukcja obsługi
Návod k použití
Kullanma kılavuzu
Használati utasítás
Manual d'instruccions
Οδηγίες χρήσης
شاميل عتال ليلد

Instrucciones de seguridad	5
Safety instructions	6
Instructions de sécurité	8
Sicherheitshinweise	9
Istruzioni di sicurezza	11
Instruções de segurança	13
Veiligheidsinstructies	14
Instrukcje bezpieczeństwa	16
Bezpečnostní pokyny	18
Güvenlik talimatları	19
Biztonsági utasítások	20
Instruccions de seguretat	22
Οδηγίες ασφαλείας	24
تعليمات السلامة	25

ÍNDICE

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SAFETY INSTRUCTIONS

- This appliance can be used by children aged 8 years and above and people with reduced physical, sensory, or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning the use of the appliance in a safe way and understand the hazards involved. Children must not play with the appliance. Cleaning and user maintenance must not be carried out by children without supervision.
- This appliance is designed for domestic use only and is not

intended for use in bars, restaurants, farmhouses, hotels, motels and offices.

- The device must only be supplied at a safety extra-low voltage in accordance with the power indicated on the device
- The device must be supplied via a residual current device (RCD) with a rated residual operating current not exceeding 30 mA.
- This appliance is rated for use in countries with warm and uniformly humid climates. It is also suitable for use in other countries.
- This device has irreplaceable batteries.

Battery instructions

- This appliance contains an irreplaceable battery.
- This appliance includes a Li-ion battery. Do not incinerate or expose it to high temperatures, as it may explode.
- Leaks from the battery or cells can occur under extreme conditions. Do not touch any liquid that leaks from the battery. If the liquid comes into contact with skin, wash immediately with soap and water. If the liquid gets into the eyes, wash them immediately with clean water for a minimum of 10 minutes and seek medical attention. Wear gloves to handle the battery and dispose of it immediately in accordance with local regulations.
- Avoid contact between the battery and small metallic objects such as paper clips, coins, keys, nails, or screws.
- Before a new battery is used or the battery is used for the first time after long-term storage, fully charge the battery.
- When the device is stored and not used for long periods of time, ensure that the battery pack is fully charged.
- Keep batteries out of the reach of children when there is a risk of ingestion.

- Battery ingestion can cause burns, soft tissue perforation and death. Severe burns may occur within two hours of ingestion.

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1. PARTS AND COMPONENTS

Fig. 1

1. Intensity level
2. Battery level
3. Intensity decrease button
4. Power button (On/Off)
5. Charging port
6. Intensity decrease button
7. Ball head attachment
8. Y-shaped head attachment
9. Thumb head attachment
10. Spine head attachment
11. Small flat head attachment
12. Conical head attachment

NOTE:

The graphics in this manual are schematic representations and may not exactly match the product

2. BEFORE USE

- Take the appliance out of the box.
- Remove all packaging materials. Keep the original box.
-
- This appliance comes in a packaging designed to protect it during transport. Take the appliance out of its box and remove all packaging materials. You can keep the original box and other packaging materials in a safe place to prevent damage to the appliance if you need to transport it in the future. If you wish to dispose of the original packaging, make sure all items are recycled properly.
- Check that all parts and components are included and in good condition. If any of them are missing or damaged, please contact Cecotec's Official Technical Support Service immediately.

Box content:

- SkinCare Massage Gun
- Charging cable
- Ball head attachment
- Y-shaped head attachment
- Spine head attachment

ENGLISH

- Small flat head attachment
 - Thumb head attachment
 - Conical head attachment
 - Instruction manual
-
- Before first use, fully discharge and then fully charge the device.
 - Do not remove the product's serial number in order to keep proper traceability if technical assistance is required.

3. OPERATION

- Insert the desired attachment or head into the massage gun's socket, ensuring it is securely fitted.
- To remove the attachment, pull it straight out.
- Avoid massaging directly on bones.

Attachment usage guide

Ball head:

Suitable for massaging large muscle groups such as arms and back.

Y-shaped head:

Suitable for massaging shoulders and areas alongside the spine to simulate manual massage.

Spine head:

Suitable for targeting deep meridian points, joints, palms, soles, etc.

Small flat head:

Suitable for full-body muscle relaxation massage.

Thumb head:

Suitable for massaging the thenar muscles of the palm and finger tendons.

Conical head:

Suitable for body meridian points.

Switching On/Off

To turn the massage gun on or off, press and hold the power button for a few seconds.

Increasing/decreasing intensity

- To increase the intensity, press/hold down the intensity increase or + button. If you hold down the button, the intensity will increase faster.
- To decrease the intensity, press/hold down the intensity decrease or - button. If you hold down the button, the intensity will decrease faster.

4. CLEANING AND MAINTENANCE

- Switch off and unplug the device from the mains before cleaning or assembling/ disassembling any parts.
- Do not immerse the device body in water or any other liquid.
- The accessories are water-resistant; therefore, they can be washed with water if necessary. Dry them thoroughly before storing or using them again.

5. DISPOSAL OF OLD ELECTRICAL AND ELECTRONIC APPLIANCES

This symbol indicates that, according to the applicable regulations, the product and/or batteries must be disposed of separately from household waste. When this product reaches the end of its service life, you must remove the batteries/ accumulators and take it to a collection point designated by local authorities.

For detailed information on how to properly dispose of electrical and electronic equipment and/or the corresponding batteries, consumers should contact their local authorities.

Compliance with the above guidelines will help to protect the environment.

6. COPYRIGHT

The intellectual property rights over the texts in this manual belong to CECOTEC INNOVACIONES, S.L. All rights reserved. The content of this publication may not, either in part or in its entirety, be reproduced, stored in a retrieval system, transmitted or distributed by any means (electronic, mechanical, photocopying, recording or similar) without prior authorisation from CECOTEC INNOVACIONES, S.L.

7. SIMPLIFIED EU DECLARATION OF CONFORMITY

Cecotec Innovaciones hereby declares that this product complies with the essential requirements and other relevant provisions of the regulations applicable in the European Union. This appliance has been designed, manufactured and tested in compliance with required safety and quality standards. The full text of the EU Declaration of Conformity can be found on the following website: <https://cecotec.es/es/information/declaration-of-conformity>

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