cecotec

SURFACE PRECISION 9500 SMARTH HEALTHY

SURFACE PRECISION 9600 SMARTH HEALTHY

Báscula de baño / Bathroom scale

Manual de instrucciones
Instruction manual
Manuel d'instructions
Bedienungsanleitung
Manuale di istruzioni
Manual de instruções
Instrukcja obsługi
Návod k použití

Instrucciones de seguridad	3
Safety instructions	5
Instructions de sécurité	8
Sicherheitshinweise	10
Istruzioni di sicurezza	13
Instruções de segurança	16
Instrukcje bezpieczeństwa	18
Bezpečnostní pokyny	21

ÍNDICE

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SAFETY INSTRUCTIONS

Read these instructions thoroughly before using the appliance. Keep this instruction manual for future reference or new users. This device is designed for personal use only and is not suitable for commercial use.

Do not use the device if it shows visible signs of damage. In case the device does not operate correctly, please contact our Technical Support Service.

Keep the device away from electromagnetic fields such as mobile telephones, remote controls, etc. Measurement accuracy could be affected.

Keep the device away from heat sources.

Do not operate the device outdoors.

Do not expose the device to extreme temperatures, moisture, wet conditions or direct sunlight.

Keep water and other liquids away from the device.

Wide temperature variations could cause condensation. In this case, allow the device to go back to room temperature before use.

Do not drop the device or jump on its measuring surface.

Do not stand on the scale with wet feet and make sure the measuring surface and the surface underneath is dry.

Do not stand on the edge of the scale.

Do not use caustic or abrasive cleaning materials to clean the device, this could damage the surface.

This scale is intended to be used for determining bodyweight, analysing water, body fat, muscle and bone mass proportion in the body, body age, visceral fat, metabolism, physical conditions and calculating Body Mass Index (BMI).

BMI function is only suitable for calculating the ideal weight of adults aged 18 and above.

The calculated results are simple reference values. Do not use these results as the basis for any form of medical treatment or diet and contact a professional for any doubt.

Body analysis function is not suitable for:

People under the age of 10 or over the age of 70.

People with a pacemaker or any other medical implant with electronic components.

People with fever, oedema or osteoporosis symptoms.

People undergoing dialysis treatment.

People who are taking cardio-vascular medication.

Pregnant women.

The device withstands a maximum body weight of 180 kg.

Keep the device, as well as its packaging material, out of the

reach of children and unsupervised animals.

Battery instructions

Keep batteries away from children and animals, they could be deadly if swallowed. In case of ingestion, seek medical assistance immediately.

Should a battery have leaked, do not allow your skin, eyes or mucous membranes to come into contact with the battery acid. In case of contact with acid, rinse the affected areas with abundant fresh water immediately and contact a doctor as soon as possible.

Make sure the battery contacts in the device are clean before inserting batteries; if they are not, clean them.

Use only the battery type specified in the technical specifications. Always replace all the batteries at the same time. Do not use different types of batteries or of different capacities.

When replacing batteries, note the correct polarity (+/-).

Take the batteries out of the device when they are discharged or when the device is not going to be used in a long time.

Batteries must not be recharged or reactivated, dismantled, thrown into fire, immersed in fluids or short-circuited.

This device requires 4 x 1.5-volt batteries, AAA type.

Warning: If "Lo" appears on the display, batteries need to be replaced.

Replacing batteries:

Open battery compartment on the bottom part of the scale.

Remove used batteries.

Introduce 4 x 1.5-volt batteries, AAA type. Pay attention to polarity.

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1. PARTS AND COMPONENTS

Fig. 1

1. Display
2. Right foot contact area
3. Left foot contact area
4. High-precision sensor
5. Battery cover
6. Body measuring tape

2. BEFORE USE

Take the product out of the box. Check carefully if there is any damage, and if there is, immediately notify Cecotec's Technical Support Service.

Remove the protective strip from the battery compartment, located on the bottom part of the device, in order to activate it. This protective strip is placed for a more secure transportation.

Further information

Weight yourself without any clothing and barefoot.

It is suggested to weigh yourself at the same time and always under the same circumstances. Measured values might be distorted if weight is measured immediately after doing sport or if you have not eaten or drunk enough over the day.

Measurement of body fat, water content and muscle and bone mass is only possible with bare and dry feet.

The scale reacts to a minimum weight of 5 kg and rounds up weight to whole 100 g.

Measuring surface must be clean in order to work properly.

3. OPERATION

Normal mode

Place the device on a solid, dry and flat surface (avoid carpets). Uneven floors will result in inaccurate measurement results.

Tap your foot on the scale in order to activate the display.

Wait until display shows '0.0' (kg/lb/st).

Before standing on the scale, select unit of measure: Kg, lb or st. See the Measuring unit setting section.

ENGLISH

Stand on the scale, the device will start calculating weight.

Once the weight has been established, the number on the display will flash. You can now get off the scale.

The device will automatically turn off after a few seconds.

Babyweight function*

This function is designed to automatically obtain precise small weight results.

This function allows seeing, activating and checking the evolution of weight from the App (see mobile APP section).

*This option is also suitable for weighing pets and baggages.

4. MOBILE APP

Scan the QR code to download and install the last version of the Surface Precision App, or browse it on Google App or App Store. Open the app and follow the instructions to link and use the scale.

Fig. 12

Register

Register in the App by pressing on "Register". On the next display, introduce your mail and a password.

Fig. 2

Log in

To log in, enter your mail and the registered account.

Fig. 3

Create profile

Enter the profile information, you can also choose the athlete mode for athletes.

Fig. 4

Linking

Link the scale via Bluetooth. Click on "Account", "Device" and "+" to search and link the scale.

Once successfully connected, the Bluetooth icon will appear on the scale and the word "Connected" will appear at the top right of the main application screen.

Fig. 5

Main display

From this screen you can see your current weight, the target weight and also, see more details about the weight.

Fig. 6

Baby mode

From the main screen you can access the baby mode, this mode will allow you to weigh your baby. Follow the instructions on the App in order to do this.

Fig. 7

Graphics

From the graphics section you can see your weights' record, which you can filter by:Recent weights

Weekly

Monthly

Yearly

You can also compare weights by clicking on the clock at the top right cornerFig. 8

Weight comparison

From this section you can compare your weight. To do so, click on the "+" on the top right corner and choose the two days you wish to compare.

Fig. 9

Settings

From this section you can change the password, manage the users, access the App settings and access the help section.

Fig. 10

Fig. 11

5. WARNING NOTIFICATIONS

The following warning notifications may appear on the display:

Display	Meaning
Lo	The batteries in the scale are too low. Insert new batteries (see the "Replacing the batteries" section).
Err	The weigh on the scale is over 180 kg, the scale does not withstand the weight.
C	An error has occurred during measurement. Step off the scale and step on it again to repeat the measuring process.

Possible errors

Error	Possible meaning	
Non-plausible values for body fat and water content.	Body fat includes relatively little water. This means the body's water content may be below the guideline figures in people who have a high body fat content. In case of endurance or performance athletes, by contrast, the guideline figures may be exceeded due to the lower fat contents and the high proportion of muscle.	
Incorrect weight measurement.	The scale is positioned on an uneven surface (e.g. a carpet). The measured values may be distorted if the surface underneath is uneven.	
	No values for body fat and water content can be calculated. Feet must be properly placed on the feet contact area, bare and dry.	

6. CLEANING AND MAINTENANCE

Do not immerse the whole unit into water or other liquids. Never immerse the device in water or other liquids.

Do not use any caustic or abrasive cleaning agents to clean the device.

Clean the device with a slightly dampened, soft cloth and, if necessary, a mild cleaning agent.

When not in use, store the device in a dry, clean place which is protected from direct sunlight and is out of the reach of children.

7. TECHNICAL SPECIFICATIONS

04090 Surface Precision 9500 Smart Healthy/04091 Surface Precision 9600 Smart Healthy

Maximum load: 180 kg

Supply: 6 V (4xAAA batteries)

Size:

9500: 310 x 310 x 20 mm

9600: 310 x 310 x 20 mm

Measuring units: Kilo (kg) / Stone (st) / Libra (lb) Net weight:

9500: 2 kg 9600: 1,86 kg Display: 82*44 mm Users: 8 Functions: 10

8. DISPOSAL OF OLD ELECTRICAL APPLIANCES

The European directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE), specifies that old household electrical appliances must not be disposed of with the normal unsorted municipal waste. Old appliances must be collected separately, in order to optimise the recovery and recycling of the materials they contain, and reduce the impact on human health and the environment.

The crossed out "wheeled bin" symbol on the product reminds you of your obligation to dispose of the appliance correctly.

Consumers must contact their local authorities or retailer for information concerning the correct disposal of old appliances and/or their batteries.

9. TECHNICAL SUPPORT AND WARRANTY

This product is under warranty for 2 years from the date of purchase, as long as the proof of purchase is submitted, the product is in perfect physical condition, and it has been given proper use, as explained in this instruction manual.

The warranty will not cover the following situations:

The product has been used for purposes other than those intended for it, misused, beaten, exposed to moisture, immersed in liquid or corrosive substances, as well as any other fault attributable to the customer.

The product has been disassembled, modified, or repaired by persons, not authorised by the official Technical Support Service of Cecotec.

Faults deriving from the normal wear and tear of its parts, due to use.

The warranty service covers every manufacturing defects of your appliance for 2 years, based on current legislation, except consumable parts. In the event of misuse, the warranty will not apply.

If at any moment you detect any problem with your product or have any doubt, do not hesitate to contact Cecotec Technical Support Service at +34 963 210 728.

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